


bkViewer

The bkViewer¹ is a stand-alone viewer for 3D ultrasound data acquired using BK Medical ultrasound systems

- Easy-to-use interface based on BK Medical ultrasound systems
- Compatible with old and new versions of BK Medical ultrasound systems' software
- Save 3D volumes in DICOM compliant format
- Retrieve and store volumes on compatible PACS systems
- Multislice view for an easy way to step through a 3D volume



bkViewer program running on a Mac

¹The bkViewer is not for diagnostic use. It is not classified as a medical device and should only be used for training and education.

Specifications

Product Name	bkViewer
Reference Type	7000
License	UA7001
Compatible Ultrasound Systems	bkSpecto ¹ , bk3000 ² , bk5000 ²
Application Software	Desktop computer ³
3D Acquisition Type	Motorized, Freehand ⁴
Display Modes	3D Cube, Render, MIP (Maximum intensity projection), 4UP, Multislice
Imaging Modes	B-mode, Color mode
3D Review	Measurements, Annotations, Volume rotation, Plane sliding & angling, Hue color, Reposition, Contrast & Brightness
3D Features	Image capture, Auto capture, Anonymize , Orientation, Wireframe, Stitch line ⁵
Supported File Types	L3D ⁶ , DICOM (with Encrypted DICOM Transfers)

¹ bkViewer is only compatible with application packs UA2639 Colorectal and UA2640 Pelvic floor.

² bkViewer is only compatible with application packs UA2644 Colorectal and UA2645 Pelvic floor.

³ See operating system requirements.

⁴ A 3D data set acquired using freehand scans cannot be used to make accurate measurements.

⁵ Only for volumes acquired with a motorized fan acquisition X14L4 (9038) transducer.

⁶ BK Medical proprietary format.

<p>SAFETY Designed and tested in accordance with EN 62304 (IEC 62304): 2006 “Medical Device Software: Software Life-cycle Processes”</p> <p>OPERATING SYSTEM REQUIREMENTS Windows: 10 or later MacOS: 10.12 or later</p> <p>HARDWARE REQUIREMENTS CPU: Intel i5 5th generation or better GPU: Intel HD Graphics 5500 or better Display: 1366x 768 pixels Memory: 8GB memory or more</p>	<p>OTHER FEATURES The about page in the program holds additional general information.</p> <p>DICOM Query/Retrieve Reports 3D Transparency¹ 3D Threshold¹ Color Transparency¹ Color MIP²</p> <p>¹ Only operational in render mode. ² Only operational in MIP mode.</p>	
---	--	--

Should any serious incident occur in relation to this device, you should report this to the manufacturer.

English source version
16-500122-EN-00

 BK Medical ApS, Mileparken 34, 2730 Herlev, Denmark. T +45 4452 8100 F +45 4452 8199

bk medical 
a GE Healthcare company

**North America
Sales and Service**
BK Medical
25 Corporate Drive,
2nd Floor, Suite 230
Burlington, MA 01803
USA
T + 1 978-326-1300
bkmedical.com

**Europe and Rest of World
Sales, Service & Design Center**
BK Medical
Mileparken 34
2730 Herlev
Denmark
T +45 4452 8100
bkmedical.com